Bedienungsanleitung
Instruction manual
Manual de instrucciones
Mode d'emploi

Togo


DE

EN

ES

FR

Enjoy mobility.

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Table of content

1.	Preliminary remark	22	5.	Using the pushing aid	29
1.1.	Signs and symbols	22	5.1.	Functional description of the control unit	29
2.	Important safety information	23	5.2.	Attaching the pushing aid	30
2.1.	General safety information	23	5.3.	Removing the pushing aid	31
2.2.	Advice on EMC disturbances	24	5.4.	Driving	31
2.3.	Braking	24	5.5.	Pushing mode	32
2.4.	Safety in vehicles and transport regulations	25	6.	Technical data	32
2.5.	Participation in traffic	25	7.	Advice on product care	33
2.6.	Purpose	25	7.1.	Cleaning and disinfection	33
2.7.	Indications	26	7.2.	Checks prior to use	34
2.8.	Contraindications	26	7.2.1.	Checking the brakes	34
2.9.	Declaration of conformity	26	7.2.2.	Checking the tyres	34
3.	Product and delivery overview	26	7.3.	Lithium-ion battery pack	34
3.1.	Inspection of delivery	26	7.4.	Charging the battery pack	34
3.2.	Scope of delivery	27	8.	Service notes for specialist retailers	35
3.3.	Overview	27	9.	Disposal	35
3.4.	Nameplate and serial number	28	10.	Transferring the pushing aid	36
4.	Installation	29	11.	Maintenance	36
4.1.	Push handle height	29	12.	Warranty	36

1. Preliminary remark

Dear user:

You have chosen a top quality electric pushing aid from Bischoff & Bischoff. We thank you for your trust in us.

The pushing aid follows the typically simple B+B use and operation concept. The pushing aid can be folded for transport in a car. The Togo pushing aid fits most current wheelchairs thanks to its universal adapter plates.

All functions, such as on/off, forwards/reverse, speed regulation and control, can be operated from the central operating unit on the handlebar. A wheelchair equipped with the Togo pushing aid can optionally be provided with active lighting, making it well-suited to use outdoors.

Before operating for the first time, please read this manual carefully. This manual is a permanent and necessary part of the Togo pushing aid. Please keep the instruction manual within reach and place it at the new user's disposal if the pushing aid is transferred to another user.

Repair work and adjustments require specific technical training and, therefore, may only be carried out by specialist retailers authorised by Bischoff & Bischoff.

Caution! Indicates information that is particularly relevant to safety. Follow the instructions for use.

Important! Indicates information that is particularly useful in the respective context.

2. Important safety information

2.1. General safety information

- In order to avoid toppling and dangerous situations, you should practice safe handling of your new wheelchair drive on a flat and familiar surface. You will be shown how to use the pushing aid and all functions by an authorised specialist retailer.
- Mounting the bracket onto the standard wheelchair must be undertaken by an authorised specialist retailer.
- Prior to using the pushing aid, please check whether all parts attached have been fastened properly.
- The user (pusher) must be mentally and physically capable of operating the Togo pushing aid.
- The values given in the operating instructions for the wheelchair, in particular air pressure, max. incline and max. user load must be adhered to.
- The pushing aid may only be operated by one assisting person (pusher).
- The wheelchair user (user) may not use the Togo pushing aid themselves as a wheelchair drive unit.
- Please use ramps to avoid obstacles where possible.
- Anti-tippers are a requirement when moving on inclines of 8% or more.
- The maximum incline allowed is 10%.
- Ledges or kerbs should only be negotiated with the pushing aid drive disconnected or after removing the pushing aid.
- The control unit must always be switched off when leaving or getting into the wheelchair.
- Observe the road traffic regulations when on the road.
- The optional lighting on the wheelchair must at all times be visible to other road users.
- If possible, wear bright, noticeable clothes. You will then be seen better by other traffic participants.
- After use, switch off the pushing aid and secure the Togo against unauthorised use by removing the key.

- Reporting incidents. The operator or user must immediately report
 1. any malfunction
 2. any change to the features or performance as well as
 3. any inaccuracies in the marking or operating instructions of a medical device, which has or could have led to death or to a serious deterioration of the health of a patient, an employee or a third party, to the Federal Institute for Drugs and Medical Devices (BfArM). This institute will immediately forward the notification to the operator's competent authority and also inform the manufacturer and the manufacturer's competent authority.

Please note that disconnecting the driving mechanisms to push the wheelchair may only be undertaken on a horizontal surface. As a result, the brake is put out of operation. Danger to life in case of non-observance!

2.2. Advice on EMC disturbances

Despite compliance with all the applicable EMC regulations and standards, it is possible that the electric drive may be disrupted by other electrical devices (e.g. mobile phones) or that it may disrupt such devices. If you observe such behaviour, switch off either the external device or your pushing aid and move the pushing aid out of reach of the disturbance.

2.3. Braking

Braking until a complete stop with the pushing aid active is effected by the automatic brake after having released the control lever.

No locking brake is available in pushing mode when the driving mechanisms are disconnected. The drives must be reconnected to lock the brakes.

2.4. Safety in vehicles and transport regulations

For transporting in vehicles (cars or similar), remove the pushing aid from the wheelchair and secure the wheelchair in accordance with the manufacturer's instructions (see operating instructions).

The pushing aid can be folded together for storage in the luggage space. In open transport spaces, secure the pushing aid against sliding using straps.

The battery pack must be at least 50% discharged for transporting in, for example, aircraft.

2.5. Participation in traffic

The pushing aid is designed for both indoor and outdoor use. Please note that you are participating in public traffic and must comply with road traffic regulations.

Do not endanger other participants through careless driving, particularly on footpaths.

2.6. Purpose

The Togo pushing aid is an electrical auxiliary drive for manual wheelchairs and is designed for both indoor and outdoor use.

The Togo pushing aid is thus intended solely for use with third parties, i.e. to support the assisting person (the pusher) when pushing.

The maximum permitted overall weight is the max. user load of the wheelchair plus the weight of the wheelchair, up to a limit of 180 kg in total.

2.7. Indications

The Togo pushing aid supports the assisting person when pushing and braking the wheelchair. Using the Togo pushing aid allows extended trips to be managed again. The use of the Togo pushing aid is indicated when the wheelchair user cannot move the wheelchair themselves and the assisting person does not have sufficient strength to push the wheelchair.

2.8. Contraindications

The assisting person must be physically and mentally in the position to safely operate the pushing aid in all operating situations.

2.9. Declaration of conformity

As the manufacturer with sole responsibility, Bischoff & Bischoff declares that the Togo pushing aid conforms to the specifications of directive 93/42/EEC.

3. Product and delivery overview

3.1. Inspection of delivery

All B+B products are subjected to a proper, final, in-house inspection and provided with the CE mark. The wheelchair is delivered in a special cardboard box. Please keep the cardboard box, if possible, after having unpacked the wheelchair. It may serve for storing the wheelchair at a later time or returning the product.

Supply, fitting the adapter and instructions will generally be undertaken by the qualified medical supplies retailer.

If delivered by rail or a forwarding agent, the product is to be checked with regard to transportation damage (cardboard box) at once in the presence of the shipper. If damaged, please contact our customer service immediately (phone: +49 700/ 6000 7070).

Please check the content with regard to completeness and sound condition. In the case of irregularities or damage, please contact our customer service (phone: +49 700/6000 7070).

3.2. Scope of delivery

Upon receipt, immediately check the content for completeness. The content consists of:


- Drive unit
- Battery pack + Charger
- Adapter plates (optional)
- Operating instructions for pushing aid and charger

3.3. Overview

3.4. Nameplate and serial number

The nameplate (fig. 02) and serial number are located on the rear of the recess for the battery unit.

- A** Manufacturer logo
- B** Manufacturer information
- C** RoHS-Logo
- D** Caution! - Follow instructions for us
- E** CE-mark
- F** Not suitable for use as a seat in a vehicle
- G** Important - Follow instructions for use
- H** Model designation
- I** Model number
- J** Maximum user weight
- K** Date of manufacture
- L** Maximum laden weight
- M** UDI number
- N** Serial number

4. Installation

The adapter plate may only be attached by a qualified medical supplies retailer.

4.1. Push handle height

The height of the push handle can be adapted to the requirements of the assisting person.

- Unscrew the plastic nut and pull the bolt completely out (fig. 04).
- Now you can change the height of the push handle.
- Push the bolt into the nearest hole and firmly tighten the plastic nut.

5. Using the pushing aid

5.1. Functional description of the control unit

Having adapted the pushing aid to your personal needs, make yourself familiar with the functions of the control unit. All control elements can be operated by an assisting person standing behind the pushing aid.

Key switch

The pushing aid is switched on and off by operating the key (fig. 05). In addition, the lockable control unit offers reliable protection against misuse.


Speed levels


The speed can be preset at 6 levels on the rotary switch (fig. 06). For safety reasons, there is only one level for reversing, at approx. 1.5 km/h.

Control lever

The electrical drive of the pushing aid is activated by operating the control lever (fig. 07).

Battery level indicator


The lighting of the battery indicator indicates that the unit is ready for operation and shows the current state of the batteries (fig. 08).

All LEDs light - max. charge

3 LEDs light - charging necessary

Pushing mode

The electrical drive is disconnected by pushing the lever downwards (fig. 09). In this way, the pushing aid can be pushed without the additional electrical drive.

Caution! The braking function is non-operational as a result.

After the lever is pulled back to the initial position, the drive engages again.

5.2. Attaching the pushing aid

Thanks to its quick-release fasteners, the Togo pushing aid can easily be attached to the wheelchair without using tools. Ideally, the wheelchair user remains in the wheelchair during attachment.

- Apply the locking brake on the wheelchair.
- Switch off the drive unit on the pushing aid and unlock it for pushing mode.
- Ensure that the locking mechanism is open.
- Move the pushing aid into the lower guide until it stops. The Togo must be slightly tilted for this (fig. 10).
- Now lever the pushing aid into the upper guides (fig. 11).
- Locking occurs separately on the right and left. To do so, pull the lever (to release) and then swing it upwards.
- Check that both sides are properly locked.

5.3. Removing the pushing aid

Removing the pushing aid is as simple as attaching it. The wheelchair user can also remain seated for this.

- Apply the locking brake on the wheelchair.
- Switch off the drive unit on the pushing aid and unlock it for pushing mode.
- Open the locking mechanism.
- Now lever the pushing aid out of the upper guide (fig. 12).
- You can now pull the pushing aid completely out of the bracket (fig. 13).
- The pushing aid removed can be folded for storage, for example in a car luggage compartment, using a simple-to-operate lever (fig. 14).



5.4. Driving

In order to avoid toppling and dangerous situations, you should initially practice safe handling of your new Togo pushing aid on a flat and familiar surface.

- First select the speed (please use the lowest level to practice) and the direction on the speed level switch.
- Switch on the pushing aid with the key switch.
- Release the locking brake on the wheelchair.
- The Togo starts smoothly by operating the control lever. Braking occurs automatically by releasing the control lever.
- The wheelchair is steered by setting the direction on the push handle.

This device is equipped with a surge protector. If the device fails when under strain, please turn off Togo and then turn it on again. Your Togo should then be ready for use again.

5.5. Pushing mode

For pure pushing mode, the Togo can be removed from the wheelchair as already described. Then push the wheelchair by using the push handles on the wheelchair.

You can, however, leave the drive unit on the wheelchair and disengage the drive motor.

- To do so, press the lever downwards (fig. 16).

In this way, the wheelchair can be pushed without the additional electrical drive. The wheelchair user can also move the wheelchair for short distances with the drive engaged if necessary.

6. Technical data

Measurements and weights

Max. user weight:	125 kg
Max. permitted total weight:	180 kg
Total length:	360 mm
Total width:	490 mm (push handle)
Overall height:	1070 - 1170 mm
Total weight wo. battery:	17,3 kg
Weight of the battery pack:	2,0 kg
Weight of the charger:	0,4 kg

Tyres

Drive wheels:	200 x 50 mm (solid rubber)
---------------	----------------------------

Brake

Brake:	electromagnetic
--------	-----------------

Electrics

Battery:	24 V 11.5 Ah lithium-ion battery pack
Operating voltage:	18.0 V - 29.4 V DC
Operating temperature:	-20°C to +45°C
Charger:	HP 1202B
Input:	100-240 VAC, 50/60 Hz, 1,2 A
Output:	29,4 VDC, 2,0 A

Motor:	DC 24 V / BLDC motor, RM 90 series
Lighting system:	lighting for attaching to the wheelchair in accordance with the German Road Traffic Act (optional)

Driving characteristics

Speed (km/h):	1.0 to 5.5 (6 levels) approx. 1.5 (safety limit)
Range (km)	approx. 14
Maximum safe incline (°):	6
Maximum climbing capability (°):	9

The maximum permitted overall weight is the max. user load of the wheelchair plus the weight of the wheelchair, up to a limit of 180 kg in total.

7. Advice on product care

7.1. Cleaning and disinfection

Make sure the Togo is switched off before starting to clean it.

- The frame and coverings can be wiped using a damp cloth. In the case of heavy soiling, use a mild cleaning agent.
- To clean the wheels, use a damp brush with plastic bristles (do not use a steel brush!).
- The control unit and hand grips can be disinfected with standard surface disinfectant.
- Make sure that electric or electronic components do not come into contact with water!

Caution! Do not use a high-pressure or steam cleaning device.

7.2. Checks prior to use

7.2.1. Checking the brakes

Prior to each use, please check the operational reliability of the brakes.

- Switch on the Togo, start carefully in the lowest level and then release the control lever. The pushing aid must come to a stop as usual.
- The motor brake is working correctly if you cannot push the wheelchair without operating the control lever.
- If it is nevertheless possible, the brake is defective. In such cases, you must no longer use the pushing aid. Please contact your specialist retailer at once.

7.2.2. Checking the tyres

Prior to each use, please check the sound condition of the tyres. Regularly check the profile depth of the tyres once a month. The profile depth should measure at least 1 mm. Defective or worn tyres must be replaced by a specialist retailer.

7.3. Lithium-ion battery pack

Only use charger supplied for charging.

Before first use, the battery pack must be fully charged in accordance with the following instructions.

The battery pack consists of maintenance-free batteries with no need to check the fluid state and that are protected against total discharge.

7.4. Charging the battery pack

Prior to each use, check the charge status of the batteries by looking at the battery status display of your control unit. If only the red LEDs of the battery status display are illuminated, you should no longer undertake any long trips and connect the charger to the pushing aid at

once. To prevent this situation, you can keep a second, charged battery pack available.

1. Plug the mains plug into the socket. The LED illuminates RED.
2. You can then connect the battery pack and charging begins automatically. The LED illuminates ORANGE.
3. The charging process is complete when the LED illuminates GREEN.

The battery pack can be charged while it is in the pushing aid (upper Fig.). Alternatively, the battery pack can be charged at a different location after removal from the pushing aid (lower Fig, battery pack). The charger should not remain attached to the fully charged battery pack for more than 8 hours!

8. Service notes for specialist retailers

We recommend having your pushing aid inspected by a specialist retailer once a year.

In the case of faults and defects on your pushing aid, return it to your specialist retailer immediately for repair.

On request, B+B will provide specialist retailers with required information, documents and wiring diagrams for repair and service. A training programme for specialist retailers is available.

In case of queries, please contact the B+B Customer Service at phone no. +49 700/60000 7070.

9. Disposal

If your pushing aid is not used any more and is to be disposed of, please contact your specialist retailer.

If you wish to dispose of the pushing aid yourself, please contact your local disposal companies regarding the applicable local regulations.

10. Transferring the pushing aid

When the pushing aid is transferred to another person, ensure that all technical documentation necessary for its safe use is given to the new user.

The pushing aid must be checked in accordance with the servicing schedule and be in perfect working order.

11. Maintenance

A servicing schedule can be made available to you by our customer service upon request.

12. Warranty

This warranty covers all defects in the pushing aid that are shown to be caused by material or production defects. The general warranty is 2 years, for batteries it is 6 months, for repairs it is 1 year. In the case of complaints, please return the product together with the warranty certificate completely filled in.

Caution! Non-observance of the instruction manual as well as servicing not carried out properly and in particular technical modifications and additions (fittings) without the approval of Bischoff & Bischoff result in both the warranty and the product liability in general lapsing.

Our general conditions of business apply.

Warranty certificate

(Please note our general business terms and conditions)

Model: Togo

Serial number: *

* (to be filled in by the specialist retailer)

Specialist retailer:

Date and stamp

*Information plates can be found on the side frame

Tip!

In spite of careful research and editing, incorrect information may have crept into our operating manuals. Right to technical changes reserved. Illustrations may show accessories. The current version of each of the operating manuals can be found in the download area of our website: www.bischoff-bischoff.com

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Bischoff & Bischoff GmbH
B+B France S.A.R.L.