B500 | B500-S

Packmaß/
Folding size

Elektrische
Sitzkantelung/
Electric seat tilt

Detail Felge
B500-S/
Rim B500-S

Otto Bock®

QUALITY FOR LIFE

CE

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Instructions for Use

Table of Contents

1	Introduction	84	6	Delivery and Putting into Operation	101
1.1	General Information.....	84	6.1	Delivery	101
1.2	Proper Usage.....	85	6.2	Putting into Operation	103
1.3	Liability.....	85	7	Operation.....	104
1.4	Information on Re-use	85	7.1	Adjustment Possibilities	104
1.5	Service.....	86	7.1.1	Backrest	104
2	Safety	87	7.1.2	Side Panel with Armrest	105
2.1	Standards and Directives	87	7.1.3	Control Panel	106
2.2	CE Conformity	87	7.1.4	Footrest.....	107
2.3	General Safety Instructions	87	7.2	Getting Into and Out of the Power Wheelchair.....	108
2.4	Safety Requirements for Transport, Storage and Assembly	88	7.2.1	From the Side	109
2.5	Safety Requirements for Operation	89	7.2.2	From the Front	109
2.6	Safety Requirements for Care, Maintenance and Disposal	93	7.3	Control.....	110
2.7	Requirements for the User	93	7.3.1	Control Panel	110
2.8	Safety Functions	93	7.3.2	Switching On and Off.....	112
2.9	Warning Symbols and Type Plates	94	7.3.3	Driving Function	113
3	Technical Data.....	97	7.3.4	Battery Indicator.....	114
4	Product Description.....	99	7.3.5	Electric Seat Functions.....	115
5	Transport and Storage	100	7.3.6	Drive-away Lock	117
			7.4	Releasing and Locking the Brake	117
			7.5	Batteries	119
			7.5.1	Charging.....	121
			7.5.2	Battery Charger	122
			8	Options	124
			8.1	Seat Options.....	125

8.1.1	Mechanical Seat Tilt Adjustment	126	9.2	Error.....	144
8.1.2	Electric Seat Tilt Adjustment.....	126	9.3	Defect/Failure	145
8.1.3	Mechanical Back Angle Adjustment	127	10	Maintenance, Cleaning and Care.....	147
8.1.4	Electric Back Angle Adjustment	127	10.1	Maintenance Intervals	147
8.1.5	Lift Seat Function (B500-S only).....	128	10.2	Changing the Fuse.....	150
8.1.6	Preparing Wheelchairs with Back Angle Adjustment for Transport	131	10.3	Changing Tyres	151
8.2	Lighting	132	10.4	Cleaning and Care.....	152
8.3	Suspension of Casters and Drive Wheels ..	134	11	Disposal.....	153
8.4	Front Wheel Swivel Lock.....	135			
8.5	Elevating Footrests.....	136			
8.5.1	Mechanically Elevating Footrest	136			
8.5.2	Electrically Adjustable Footrests.....	136			
8.6	Control Panel Holders	137			
8.6.1	Height Adjustable Control Panel Holder ..	137			
8.6.2	Swing-away Control Panel Holder.....	137			
8.7	External Mileage Counter	138			
8.8	Recaro® Seat.....	139			
8.9	Other Options	141			
8.9.1	Lap Belt	141			
8.9.2	Four-way Chest Strap Kit	141			
8.9.3	Headrest Mounting Kit (Standard Seat) ..	142			
8.9.4	Rear Bumper.....	142			
8.9.5	Other Options	143			
9	Failures/Troubleshooting	144			
9.1	Warning	144			

List of Figures

Fig. 1	Signage on the B500/B500-S	94	Fig. 22	Electric back angle adjustment.....	128
Fig. 2	Folding size.....	100	Fig. 23	Lift seat in raised position (B500-S).....	130
Fig. 3	Main components of the B500/B500-S... ..	102	Fig. 24	Snail symbol indicating creep speed	131
Fig. 3.1	Anti-tipper B500	103	Fig. 25	Cross bolt locked in the holding device....	131
Fig. 4	Fuse holder	103	Fig. 26	Front lights on the side panel.....	132
Fig. 5	Release strap for the backrest	104	Fig. 27	Changing the front lights	132
Fig. 6	Removing the side panel	105	Fig. 28	Replacing the direction indicator lamp.....	133
Fig. 7	Adjusting the height of the armrest	105	Fig. 29	Rear light without lamp cover	134
Fig. 8	Adapting the control panel	106	Fig. 30	Release bolt on the battery drawer.....	135
Fig. 9	Footrest lock	107	Fig. 31	Caster swivel lock engaged and released..	135
Fig. 10	Adjusting the lower leg length.....	108	Fig. 32	Mechanically elevating footrest.....	136
Fig. 11	Footplates flipped up.....	109	Fig. 33	Elevated footrest	137
Fig. 12	Footrests removed.....	109	Fig. 34	Height adjustment of control panel holder .	137
Fig. 13	Control panel	110	Fig. 35	Swing-away control panel holder.....	138
Fig. 14	Display with all its symbols	112	Fig. 36	External mileage counter	138
Fig. 15	Brake release.....	119	Fig. 37	Applying the lap belt.....	141
Fig. 16	Locking the brake.....	119	Fig. 38	Headrest mounting kit	142
Fig. 17	Snap lock on safety strap.....	120	Fig. 39	Rear bumper	143
Fig. 18	Batteries	121	Fig. 40	Fuse holder	150
Fig. 19	Mechanical seat tilt adjustment	126	Fig. 41	Disassembling drive wheel	151
Fig. 20	Electric seat tilt adjustment.....	127	Fig. 42	Disassembling front wheel.....	151
Fig. 21	Mechanical back angle adjustment	127			

1 Introduction

1.1 General Information

These instructions for use provide the user as well as his or her attendants with all the required knowledge on the design, function, operation and maintenance of the B500 power wheelchair from Otto Bock Mobility Solutions GmbH. The instructions for use contain all the information that is needed to use the power wheelchair safely, to determine the possible causes of a malfunction and to help eliminate it.

Knowledge of these instructions for use is absolutely necessary for ensuring that the power wheelchair is used safely. Therefore, the user and his or her attendants must read the instructions for use thoroughly, especially the section on "Safety", before using the power wheelchair. This will ensure the complete utilisation of all the features of the power wheelchair.

The instructions for use have been prepared on the basis of DIN EN 62079 "Preparation of instructions – Structuring, content and presentation". The instructions for use are divided into 11 sections. The header of each page contains the title of the respective section; the footer of each page contains the description of the power wheelchair, the edition of the instructions for use and the page number.

The instructions for use also contain cross references that make it easier for you to handle the manual. An example of this is "See section 4.2".

Symbols and signal words:

DANGER!

Warning of possible dangers to life and limb.

ATTENTION!

Warning of possible damage to property

NOTE!

Note on important information and tips

NOTE!

Information on environment protection

1.2 Proper Usage

The B500 power wheelchair is designed for individual mobility indoors and outdoors. It is intended solely for people who are unable to walk or who have a walking impediment. The power wheelchair may only be combined with the options mentioned in these instructions for use and vice versa. Otto Bock assumes no liability for combinations with medical products and/or accessories from other manufacturers outside of the modular system.

An option for the B500 is to have it controlled by an attendant with the use of special controls.

Using the B500 for purposes other than the ones described above is considered improper. The user, rather than the manufacturer, shall be solely liable for property damage and personal injury caused by improper use.

The B500 may only be used by persons trained to use it. Training users and attendants to use the power wheelchair is required for protecting persons from danger and for ensuring that the B500 is operated safely and correctly.

The operational safety of the B500 can only be ensured if it is used properly in accordance with the information contained in these instructions for use. The user is ultimately responsible for accident-free operation.

1.3 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

1.4 Information on Re-use

The B500 power wheelchair is suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

Based on market observations and the current state of technology, the manufacturer has calculated that the B500 power wheelchair can be used for a period of 5 years, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the wheelchair is stored at the dealer or with the cost bearer are not included in this period. It

should be clearly pointed out, however, that the B500 is a reliable product far beyond this defined period of time, provided that it is cared for and maintained appropriately.

If the product is going to be re-used, it must be thoroughly cleaned and disinfected.

Afterwards, the condition of the product must be examined by an authorised dealer to check for wear and tear as well as any damage.

Any worn and damaged components as well as components which do not fit or are unsuitable for the new user must be replaced.

The service manual includes a service plan for each model, detailed information and a list of the required tools.

1.5 Service

NOTE!

Service and repairs on the power wheelchair must only be carried out by specialist dealers with personnel authorised and trained by Otto Bock. Should any problems arise, please contact your power wheelchair supplier.

Should you have any questions or a problem that cannot be resolved despite using the instructions for use, please contact Otto Bock Customer Service (see inside cover for address).

Otto Bock endeavours to support customers in all respects in order to keep them satisfied with their product for a long time.

2 Safety

2.1 Standards and Directives

All information on safety contained in these instructions for use refers to the currently valid national laws and regulations of the European Union. In other countries compliance with the applicable laws and national regulations is required.

Besides the safety instructions contained in these instructions for use, the user is required to observe and comply with the generally valid professional association regulations, the accident prevention regulations and environment protection regulations. All information contained in these instructions for use must be complied with at all times without restrictions.

The B500 power wheelchair has been constructed in accordance with the currently valid technical rules and is safe to operate. The safety of the B500 power wheelchair has been confirmed by CE certification and the declaration of conformity.

2.2 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

2.3 General Safety Instructions

- The B500 may only be used properly.
- The B500 may only be used by persons trained in its use and must not be used by unauthorised persons.
- The B500 may be used for transporting only one person.
- For safety reasons, it is strictly recommended that the user always apply and fasten the lap belt, which is available from Otto Bock as an accessory.
- All safety instructions contained in these instructions for use and in all further applicable documents are to be observed and complied with. The instructions for use must be available to the user at all times.
- The back upholstery and seat cushion of the power wheelchair are not highly flammable, but can however catch fire. For this reason, use utmost caution near flammable objects or fire, especially lit cigarettes.

- ❑ When driving on inclines or slopes, it is not permissible to negotiate obstacles.
- ❑ Operating the power wheelchair on stairs is not permitted.
- ❑ Avoid getting into and out of the power wheelchair while standing on an incline or slope.
- ❑ If the driving mode is changed while driving, the power wheelchair will accelerate or decelerate.
- ❑ Releasing the brake may result in uncontrolled rolling of the power wheelchair. Therefore make sure that the brake is engaged after parking the power wheelchair.
- ❑ The lap belt helps to additionally stabilise the person sitting in the power wheelchair. The lap belt must never be used as part of a restraint system for transportation of the person in motor vehicles for the disabled.
- ❑ Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunction. In this case, please contact your authorised dealer immediately.

2.4 Safety Requirements for Transport, Storage and Assembly

- ❑ Only sufficiently strong hoisting devices may be used for transport.
- ❑ The power wheelchair must be secured in accordance with the requirements for the transport aid. The tensioning straps may only be attached to the corresponding fixing eyelets.
- ❑ During transport on lifting platforms, in lifts, busses, or trains the wheelchair controls must be switched off. The brake must be locked and the power wheelchair must stand in the middle of the lifting platform. None of the power wheelchair's components, such as its anti-tippers or other components, may be in the danger area.
- ❑ When the power wheelchair is not being used for an extended period of time or during shipping, remove the fuse.
- ❑ After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened.
- ❑ We recommend that, wherever and whenever possible, users of a wheelchair/ mobility base with seating shell or buggy transfer to the seats installed in the motor vehicle and use the corresponding vehicle restraint systems, because this is the only way to ensure opti-

mal protection of the passengers in case of an accident. If your B500 power wheelchair is to be used as a seating transportation in a motor vehicle the anchor brackets offered by Otto Bock and appropriate restraint systems for the wheelchair must be used. For more information, please refer to our brochure “Using your Wheelchair/ Mobility Base with Seating Shell or Buggy/Stroller for Transportation in Motor Vehicles”, order no. 646D158.

- ❑ When mounting a Recaro® seat, ensure that the seat adapter is securely attached to the retainer bushings and that the anterior locking bolts have snapped into place up to the key ring!
- ❑ When mounting a special control, the user is obliged to have a specialist dealer give him or her detailed instructions on how to use it.
- ❑ The mounting position of the seat must not be changed. Otherwise there is the risk of tipping.

2.5 Safety Requirements for Operation

- ❑ The user is obliged to make sure that the power wheelchair and its safety functions are in safe and proper condition every time before using it. The power wheelchair may only be used if all safety functions, such as the automatic brakes, are functional. Inoperable brakes can lead to severe accidents with fatal injuries.

- ❑ The power wheelchair must be inspected for functional reliability and driving safety at least once a year by an authorised specialist.
- ❑ Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents. Always test the driving characteristics of the power wheelchair after configuration / programming is complete.
- ❑ Programming must only be completed by authorised personnel. Neither Otto Bock nor the control unit manufacturer are liable for damages (especially in combination with special controls) caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.
- ❑ To ensure safe driving operation, the anti-tipper must have been mounted correctly and must be in proper condition.
- ❑ If any faults, defects or other dangers that can lead to personal injury are detected, the power wheelchair must be put out of operation immediately.
- ❑ Prior to using the power wheelchair, all the necessary mechanical adaptations (e. g. mounting special controls) and software settings (e. g. programming the

control) must be made to comply with the individual requirements and abilities of the user. The settings may only be made by trained specialists who have been authorised by Otto Bock.

- ❑ The B500 may be used within a temperature range of -25 °C to +50 °C. During operation, the wheelchair may not be exposed to temperatures below or above this range.
- ❑ The maximum load for the B500 is 140 kg/309 lbs. This load must not be exceeded.
- ❑ The lift seat function may only be used on horizontal surfaces and only with a vertical back angle.
- ❑ The seat function actuators are not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time). General recommendation: Under maximum load, a pause of approximately 90 seconds is required after 10 seconds of operation. The electric seat functions are considered independently of the driving function for this purpose.
- ❑ The lift seat may only be subjected to a maximum load of 130 kg.
- ❑ The lift seat may not be operated in case of errors or malfunctions. If creep speed is not activated when the lift seat is operated, please visit a specialist dealer immediately. The power wheelchair may only be used with the seat in its lowest position until the problem is

corrected.

- ❑ The power wheelchair must be visually checked every time before it is used to make sure the tread depth is sufficient and that the tyres are inflated to the correct tyre pressure (printed on the sidewall of the tyre).
- ❑ When the wheelchair is used on the street, its user is required to obey the public traffic regulations.
- ❑ Using the power wheelchair on very slippery ground (such as icy surfaces) is not permitted. Driving on very coarse-grained surfaces (gravel or pebbles) is not permitted either.
- ❑ When getting into or out of the power wheelchair, the controls must be switched off. The footrests and armrests are not capable of bearing full body weight, and therefore must not be used for getting into or out of the wheelchair.
- ❑ The critical obstacle height of the B500 is 5 cm. Differences in height greater than 5 cm must not be crossed.
- ❑ Risk of injury if the power wheelchair tips over! Before driving up slopes or over obstacles, a tilted seat must be lowered to its basic position and the back must be brought into the upright position.
- ❑ Drive over obstacles such as steps or curbs at reduced speed only (max. 3 km/h). Always approach obstacles at a right angle and cross over them without stopping.

- ❑ Manoeuvring is only allowed at a reduced speed.
- ❑ To lift the wheelchair, grasp it by frame components only! Never grasp the footrests or armrests to lift the wheelchair.
- ❑ All attendants must be informed that there are pinch points caused by the construction in the area between the seat frame and wheelchair frame. Users and attendants must not reach into the danger area.
- ❑ The controls of the power wheelchair are protected according to protective system IP 54 and can thus be used in bad weather (e. g. rain). The controls have been approved for indoor and outdoor use and comply with the requirements for climate and splashing water.
- ❑ The power wheelchair has been tested according to EMC regulations. The following particularities must be observed when using the wheelchair:
 - The driving characteristics of the power wheelchair can be influenced by electromagnetic fields (mobile phones or other radiating devices). For this reason, all mobile devices must be switched off when driving.
 - The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Therefore, switch off the controls whenever you do not need them.
- ❑ The tyres of the power wheelchair contain chemical substances that can react to other chemical substances (such as cleaning agents, acids, etc.).
- ❑ Direct sunlight (UV light) causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread.
- ❑ Whenever possible, avoid parking the wheelchair outdoors. Regardless of wear and tear, the tyres should be replaced every 2 years.
- ❑ Black tyres contain sooty particles, which can leave traces on the ground. For this reason, a suitable mat should be placed underneath the wheelchair, if it is not going to be used over an extended period of time.
- ❑ The B500 has been approved for driving on inclines and slopes of no more than 17 %. It is not permitted to drive on inclines or slopes steeper than this.
- ❑ As a precaution, the driving speed must be reduced when driving downhill (e.g. set it to speed level 1).
- ❑ When the brake is released, there is not any brake function available. This must be taken into account, especially when driving on slopes or inclines.
- ❑ After each emergency stop, the controls of the power wheelchair must be turned on again. In the event of communication problems in the bus system of the controls, the system triggers an emergency stop and thus prevents any uncontrolled functions. If the driving function is still not available after switching on the

controls again, the brakes must be released to switch over to the pushing mode. In this case, contact a specialist dealer as soon as possible.

- ❑ If the controls do not emit an error signal as soon as the brake lever is activated, something has malfunctioned. The settings must be checked by a specialist dealer.
- ❑ The power wheelchair can be used for transportation of the person in vehicles for the disabled, with some restrictions. Please read the information contained in section 5 (Storage/Transport).
- ❑ Explosive gases can develop while the batteries are charging. For this reason, the following precautions are to be taken during charging:
 - The controls must be turned off.
 - Provide for sufficient ventilation when using the battery charger in closed rooms.
 - Smoking and fire are not permitted. Sparks must be avoided.
 - The vent slots of the housing must not be covered.
- ❑ For power wheelchairs featuring a battery drawer: When closing the battery drawer, the release bolt must engage in the recess provided for this purpose. The release bolt must be easy to see and touch on the underside of the drive unit sustainer.

- ❑ In case of power wheelchairs which feature gas compression springs, attention must be paid to oil leaking from the piston rod. Leakages reduce the functional capabilities of the gas compression springs and result in their failure.
- ❑ Defective gas compression springs (e. g. scratches or damage to the piston rod caused by impact; deformation of the cylinder tube) must be replaced immediately. Piston rods do not need lubrication. They are maintenance-free.
- ❑ When the wheelchair is not used for extended periods of time, or if the tyres are heated strongly (e. g. near radiators or by sunlight shining through a window), permanent deformation of the tyres will result. Therefore always make sure that the wheelchair has sufficient distance from sources of heat, move your wheelchair from time to time, or jack up the wheelchair when storing it.

2.6 Safety Requirements for Care, Maintenance and Disposal

- ❑ Only trained specialists who have been authorised by Otto Bock may do maintenance work on the power wheelchair. This also applies to all repairs and settings on the brakes. Incorrect settings can lead to brake failure.
- ❑ The fuse must always be removed for any maintenance work where the battery cover is open.
- ❑ Never use a water jet or high-pressure cleaning apparatus for cleaning the power wheelchair. The water may not come into direct contact with the electronics, motor or battery.
- ❑ Defective batteries must be disposed of properly in accordance with national regulations.

2.7 Requirements for the User

- ❑ The power wheelchair may only be used by informed users. For this purpose, the user and maybe any attendant must receive instructions on how to use the power wheelchair from trained specialists who have been authorised by Otto Bock.
- ❑ The user must have read and understood all the information in the instructions for use.
- ❑ The wheelchair may not be operated in cases of

exhaustion or under the influence of alcohol or drugs. The user may not have any mental restrictions which can temporarily or permanently restrict attentiveness and judgement.

- ❑ Before using the power wheelchair for the first time, the user must become familiar with the wheelchair on level ground first. This can help to avoid potentially dangerous situations such as tipping. Getting to know how the power wheelchair reacts when the centre of gravity shifts, for example, on slopes or inclines or when clearing obstacles like steps and curbs should be done only with assistance from another person before using the wheelchair for the first time.
- ❑ When driving in the dark, it is strongly recommended that the user wear bright clothes or clothing with reflectors.

2.8 Safety Functions

In dangerous situations, the B500 can be turned off at any time using the on/off button. When the button is pressed, the power wheelchair brakes immediately.

If any malfunctions occur such as an insufficient supply of energy to the brake, the software will recognise them and trigger the emergency stop brake or reduce the speed of the power wheelchair. At the same time, an acoustic warning signal will be emitted.

2.9 Warning Symbols and Type Plates


Fig. 1 Signage on the B500/B500-S

Label/Type Plate	Explanation
<p>The image shows two labels for Otto Bock wheelchairs. The top label is for the B500 model and the bottom label is for the B500-S model. Both labels include the following information:</p> <ul style="list-style-type: none"> Brand name: Otto Bock Quality for life slogan Model name: B500 (top) or B500-S (bottom) Technical specifications: <ul style="list-style-type: none"> max. Zuladung: xxxkg (B) max. Steigfähigkeit: xx % (C) max. Geschwindigkeit: xx km/h (D) zul. Achslast vorn: xxx kg (E) zul. Achslast hinten: xxx kg (F) zul. Gesamtgewicht: xxx kg (F) Serien-Nr.: I xxxxxxxxxxxxxxxx (I) Manufacturer information: Otto Bock Mobility Solutions GmbH, Lindenstr. 13 - Königsee/Germany - www.ottoBock.com (K) Made in Germany (L) Safety symbols: G (recycling), CE (H), and J (instructions) Warning: Sicherheitshinweise beachten! 	<ul style="list-style-type: none"> A Type designation B Maximum load capacity (see Section “Technical Data”) C Maximum climbing ability (see Section “Technical Data”) D Maximum speed (see Section “Technical Data”) E Permissible axle load, front / rear F Permissible total weight G Symbol for separate collection of electric and electronic devices. Components of the power wheelchair and the batteries must not be disposed of like regular domestic waste. H CE marking – product safety in conformity with EC Directives I Serial number J Read the Instructions for Use prior to using the product. Observe the safety instructions in the Instructions for Use. K Manufacturer information / address L Country of manufacture

Label/Type Plate	Explanation
	<p>A Electric driving operation: lock motor brake B Manual pushing mode: unlock motor brake</p>
	<p>Risk of pinching. Do not reach into the danger area.</p>

Table 1 Signage on the B500 / B500-S

3 Technical Data

Measurements and weights	
Seat width	38 – 42 cm or 43 – 48 cm
Seat depth	38 – 46 cm or 42 – 50 cm
Seat height	45 – 60 cm
Armrest height	24 – 36 cm
Armrest length	26 cm
Lower leg length	25 – 34 cm or 35 – 44 cm
Back height	45 or 55 cm
Back angle	-9/1/11/21° or 0/10/20/30°
Overall width	64.5 cm
Overall height	103 cm
Overall length	108 cm
Weight empty	95 kg
Maximum load capacity	140 kg (patient weight)
Turning radius	155 cm
Tyre size	
Caster	9" – 10"
Drive wheel	14"
Air pressure	printed on the sidewall of the tyre

Electrical Installation	
Operating voltage	24 V
Batteries	
Wet cell batteries	2 x 12 V; 60 Ah (C5)
Gel batteries	2 x 12 V; 63 Ah (C5)
Control unit:	
Model	EnAble 40 with controller and control panel
Operating voltage	24 V DC
Max. output current per motor	75 A
Lighting:	
Front direction indicator	
Front light	H21W 12 V BAY9s
Rear direction indicator	HMP 08 2, 4W, 6V
Rear light	C21W; 12 V; SV8.5 C10W; 6 V; SV8.5
Fuse	80 A

Driving data	
Speed	6 km/h or 10 km/h
Climbing ability	17 %
Obstacle height that can be cleared	5 cm
Distance range	approx. 35 km
Operating temperature range	-25 °C to +50 °C
Transport and storage temperature range	-40 °C to +65 °C
Battery charger	
Model	MEG 2409 Automatic battery charger with computer-controlled characteristic curve
Power requirements	230 V -10 %; +6 %
Mains frequency	50 Hz ±4 %
Protective system	1 (protective conductor)
Charging connection	24 V DC
Nominal charging current	9 A

Battery charger	
Residual ripple	>1 %
Characteristic curve	IUU with trickle charging phase, computer-controlled
Primary fuse	T4G fuse link, not accessible from outside
Secondary fuse	Electronic, reversible reverse battery protection, short circuit-proof, idling-proof, overheating protection
Protective system	IP 21
Ambient temperature range	-10 °C to +40 °C
Indicator	2 LEDs
Weight	1.7 kg
Dimensions (W x H x D)	105 x 65 x 205 mm

Table 2 Technical data

4 Product Description

The B500 power wheelchair is suitable for indoor and outdoor use. It is compactly designed and easy to manoeuvre for use indoors. Two 12 V batteries power its high-performance drive system, which makes it possible to overcome obstacles (category B of EN 12184) and to ensure safe operation.

The control unit comprises a control panel for entering the driving commands and for displaying the current state of the functions as well as a controller that uses the input data to control the drive motors and other electric functions. Data transmission is realised via a bus system.

Since the system can be programmed, it is possible to adapt the control unit to the personal requirements of the user; e. g. the speed, acceleration and deceleration values can all be adapted.

The special features of the power wheelchair include:

- Easy servicing due to good and clear accessibility of all component groups
- Customisation through options and special fabrication with modular components (chassis, seating system, control unit, accessories)

The modular design of the B500 makes it possible to add other modules and devices in addition to the main components (see fig. 3), e. g. electric seat adjustment, tray, etc. These options are described in section 8.

The versatility of this power wheelchair and its modular design make it suitable for patients who have walking impediments or walking disabilities due to:

- Paralyses
- Loss of limbs (lower limb amputation)
- Defective or deformed limbs
- Joint contractures or defects
- Other diseases

The B500 power wheelchair was specially designed for users who are able to move independently in such a power wheelchair.

Fitting considerations:

- Body height and weight (maximum load is 140 kg / 309 lbs)
- Physical and psychological limitations
- Age of the patient
- Home conditions
- Environment

5 Transport and Storage

DANGER!

Risk of injury as a result of improper transportation!

The seats installed in the motor vehicle along with the corresponding vehicle restraint systems offer passengers ideal protection in the event of an accident. Only in exceptional cases and only when using the safety elements and restraint systems offered by Otto Bock should the B500 power wheelchair be used as a seat in motor vehicles for the disabled. For more information, please refer to our brochure "Using your Wheelchair/ Mobility Base with Seating Shell or Buggy/Stroller for Transportation in Motor Vehicles", order no. 646D158.

DANGER!

Risk of injury as a result of not engaging the locks sufficiently!

Secure the power wheelchair during transport in another vehicle by sufficiently tightening the tensioning straps.

The frame of the power wheelchair has 4 fixing eyelets for the attachment of the tensioning straps.

Before transporting the power wheelchair, switch off the controls and lock the brake.

You can reduce the size (fig. 2) of the power wheelchair for transport by folding down the backrest and removing the side panels and footrests (see section 7.1).

Fig. 2 Folding size

Tyres contain chemical substances that can react to other chemical substances (such as cleaning agents, acids, etc.).

Black tyres also contain sooty particles.

**NOTE!**

If your power wheelchair is not moved for several days, permanent colour changes might result on the ground where the wheelchair comes into contact with it. For this reason, place a suitable mat underneath the wheelchair, if it is not going to be used over an extended period of time.

**NOTE!**

Black tyres can leave black traces on the ground under some circumstances. For this reason, we recommend choosing grey tyres when using the power wheelchair mainly indoors.

**NOTE!**

Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

The power wheelchair must be stored in a dry place.

For transport and storage, an ambient temperature range from -40 °C to +65 °C must be observed.

6 Delivery and Putting into Operation

6.1 Delivery

Delivery includes:

- Adapted power wheelchair with main components (see fig. 3)
- Battery charger
- Operating instructions
- Options (see section 8)

**NOTE!**

The options included in the delivery depend on the product configuration that has been selected for the power wheelchair.

Upon delivery by the specialist dealer the power wheelchair is ready for use. All the settings correspond to the indications on the order form or are adjusted directly on site by the dealer. The power wheelchair is adapted to the personal requirements of the individual.

The functions of the individual components can be tested by following the instructions in section 7. Possible malfunctions are described in section 9.


Fig. 3 Main components of the B500 / B500-S

- | | | | |
|-----------------------------|------------------------|-----------------|--------------------------|
| 1 Backrest | 3 Armrest (side panel) | 5 Brake release | 7 Motor with drive wheel |
| 2 Joystick and control unit | 4 Seat cushion | 6 Footrest | 8 Anti-tipper |


Fig. 3.1 Anti-tipper B500

The tool kit consists of:

- 1 set of Allen wrenches, sizes 3 mm, 4 mm, 5 mm, 6 mm
- 1 open-end wrench, size 13 mm

6.2 Putting into Operation

Prior to putting the power wheelchair into operation, the completeness (see fig. 3) and function of all its components must be checked. Before the power wheelchair can be switched on, the 80 A fuse must be inserted into the fuse holder at the rear end of the battery case (see fig. 4).

To do so, remove the fuse from the protective cover on the control panel, open the cap of the safety holder and insert the fuse.

Make sure that the fuse is pressed into the centre of the spring contacts and that it does not lie at an angle. Close the cap until you feel it snap back into place.


Fig. 4 Fuse holder

- 1 Cap, open
- 2 Fuse, inserted
- 3 Fuse holder

7 Operation

7.1 Adjustment Possibilities

DANGER!

Risk of injury as a result of not engaging the lock!

Be sure to firmly retighten all screws and nuts after making adjustments.

Various adjustments can be made to the power wheelchair.

The seat height, seat width and seat angle have been set in accordance with the indications on the customer order form and may only be changed by the specialist dealer.

The following items can be adapted by the user:

- Back angle
- Armrest height
- Armrest position
- Lower leg length

If need be, the footrests and side panels can be removed.

7.1.1 Backrest

The backrest can be set to four different angles by pulling the release strap at the base of the backrest (fig. 5).

Once the backrest is in the desired position, release the strap to engage the locks.

Fig. 5 Release strap for the backrest

7.1.2 Side Panel with Armrest

The side panels can be removed by loosening the thumb screws and lifting the side panels. If the control panel is attached to the armrest, it must be disconnected before taking out the side panel (section 7.1.3).

The thumb screws are located at the lower end of the armrest holders (see arrow in fig. 6).

Fig. 6 Removing the side panel

After remounting the side panel, do not forget to firmly retighten the thumb screws.

The height of the armrest can be changed by loosening the screw (see arrow in fig. 7).


Fig. 7 Adjusting the height of the armrest

Firmly re-tighten the screws after making adjustments.

7.1.3 Control Panel

Adapting the Control Panel to Arm Length

To adapt the control panel to the arm length of the user, you must loosen the three screws on the underside of the armrest (see fig. 8).

You can now slide the rail with the control panel backwards and forwards.

Fig. 8 Adapting the control panel

Firmly re-tighten the screws after making adjustments.

NOTE!

If the control panel rail is too long, the protruding part can simply be cut off with a saw.

Mounting and Removing the Control Panel

The standard configuration has the control panel mounted on the right side. If required, it can also be mounted on the left armrest.

Spring sheet metal connects the control panel to the rail, which is attached below the armrest.

First undo the hook and loop closures that fix the cable of the control panel to the wheelchair. Then pull the control panel toward the front and off the rail. Disassemble the rail from the armrest and reattach it on the other side. Insert the control panel onto the rail and fix the cable to the wheelchair by using the hook and loop straps.

ATTENTION!

Cable damage!

Positioning the cable incorrectly can lead to pinching and thus damage to the cable. The cable must not be attached too tightly or too loosely. Avoid bending or squeezing the cable.

7.1.4 Footrest

DANGER!

Risk of getting pinched!
 Make sure your fingers are not in the danger area when flipping the footplates up or down.

Removal

1. Remove the calf band.
2. Flip the footplate up, release the footrest lock (see fig. 9) and swing the footrest either to the inside or outside.
3. Pull the footrest up to remove it.

Reassembly

1. Reinsert the footrest from above by placing it in the holding device and swinging it forward until the footrest lock engages.
2. Reattach the calf band to the holding device.

Fig. 9 Footrest lock

Adjusting the Lower Leg Length (see fig. 10)

1. Loosen the screws on the footrest bar.
2. Move the footplate up or down to adapt the height to the individual lower leg length and seat cushion thickness.

NOTE!

The footrest bar must not be pulled out of the holding device by more than 160 mm.

3. Retighten the screws.

Fig. 10 Adjusting the lower leg length

7.2 Getting Into and Out of the Power Wheelchair

The modular design of the power wheelchair and the ease with which you can remove the side panels and footrests make it easy to get into and out of the wheelchair from the side or from the front.

Getting into and out of the wheelchair can be done by the user individually in a way that suits him or her best.

DANGER!

Risk of injury if the power wheelchair starts rolling!

When getting into or out of the power wheelchair, the wheelchair controls must always be switched off.

ATTENTION!

Risk of breaking the footplates and armrests as a result of overloading them!

The user must not apply his or her entire body weight to the footrests and armrests for support when getting into or out of the wheelchair. Do not use the footplates when getting into or out of the wheelchair.

7.2.1 From the Side

Getting into the wheelchair from the side requires the user to remove the right or left side panel depending on the side where the user will get into the wheelchair.

Bring the power wheelchair as close as possible to the place where its user is sitting.

If the control panel is on the side the user wants to use to get into/out of the wheelchair, then undo the hook and loop closures that fix the cable for the control panel to the wheelchair and pull off the control panel. Remove the side panel as described in section 7.1.2. If necessary, remove the footrest (see section 7.1.4).

The user can then slide onto the wheelchair's seat from the side. Use of a transfer board will make this even easier.

7.2.2 From the Front

Flipping up the two footplates makes it possible for the user to get into and out of the wheelchair from the front (see fig. 11). The space for getting into/out of the wheelchair can be increased by removing the footrests (fig. 12).

The assistance of an attendant or a transfer lifter make it easy for the user to get into or out of the power wheelchair. Use of a rotation plate is also possible.

Fig. 11 Footplates flipped up

Fig. 12 Footrests removed

7.3 Control

ATTENTION!

Risk of impairing the driving performance of the power wheelchair!

The driving characteristics of the power wheelchair can be influenced by electromagnetic fields (mobile phones or other radiating devices). For this reason, all mobile devices are to be switched off when driving.

ATTENTION!

Risk of damage to other devices!

The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Therefore, switch off the controls whenever you do not need them.

7.3.1 Control Panel

The control panel is used to control the power wheelchair. The control panel consists of a button section, LCD display and joystick. The charging/programming receptacle is on the underside.

The control panel is used to switch the power wheelchair on and off, to enter driving commands and to display the current state of certain functions and components.


Fig. 13 Control panel

- | | |
|-----------------------------|------------------------------------|
| 1 Joystick | 6 Mode button |
| 2 Direction indicator, left | 7 Horn |
| 3 Warning flasher | 8 Direction indicator, right |
| 4 LCD display | 9 Light |
| 5 On/Off button | 10 Charging receptacle (underside) |

On/Off Button

The on/off button is used for switching the power wheelchair on, for activating the drive-away lock and for switching the power wheelchair off.

Mode Button

Pressing the mode button shortly increases the speed level. After reaching the maximum speed level, you can change back to speed level 1 by pressing the button again. Pressing the button for at least 2 seconds switches the control unit to the “Electric seat functions” menu (see section 7.3.5).

Horn

The horn will sound as long as the horn button is being pressed.

Warning Flasher

Pressing the warning flasher button activates all four direction indicators.

Light

The front and rear lights are activated and deactivated by pressing the light button.

Direction Indicators

The right or left front and rear direction indicators are activated and deactivated by pressing the “direction indicator, right” or “direction indicator, left” buttons respectively. The direction indicators automatically go out after 20 seconds.

LCD Display

The LCD display is the communication interface between the user and the control unit. It indicates the selected speed level, the remaining battery capacity, the status of electric options and special functions as well as warnings and errors. During the start phase, all display symbols are visible.


Fig. 14 Display with all its symbols

- | | |
|-----------------------------|-------------------------------|
| 1 Direction indicator, left | 9 Excess temperature |
| 2 Battery capacity | 10 Direction indicator, right |
| 3 Electric backrest | 11 Drive wheel brake |
| 4 Electric seat tilt | 12 Warning |
| 5 Control panel | 13 Creep speed |
| 6 Drive-away lock | 14 Power module |
| 7 Speed level | 15 Drive motor |
| 8 Open end wrench | 16 Light |

7.3.2 Switching On and Off

DANGER!

Danger to life if the brake fails to function!
The brake release lever must be locked when using the power wheelchair. The automatic brakes must be ready for operation and functional.

DANGER!

Risk of accident as a result of insufficient or excessive tyre pressure!
The power wheelchair must be visually checked every time before it is used to make sure that the tread depth is sufficient and that the tyres are inflated to the correct tyre pressure. Incorrect tyre pressure reduces the tyres' service life and makes the wheelchair's driving behaviour worse.

The control unit of the power wheelchair is switched on and off by pressing the on/off button (see fig. 13, item 5). If the control is not used for some time, the wheelchair is switched off automatically. It is also possible to switch off the power wheelchair with the on/off button during driving. In this case, the wheelchair brakes immediately until it stops.

NOTE!

Every time you switch on the control unit, it will be at the speed level which was selected last.

7.3.3 Driving Function

DANGER!

Risk of injury if the power wheelchair gets out of control!

Switch off the controls on the power wheelchair whenever you do not need them. This will prevent you from unintentionally activating the joystick.

DANGER!

Risk of injury if the power wheelchair tips over! Observe the following safety instructions when driving the power wheelchair:

- Do not drive on inclines and slopes of more than 17 %.
- When driving downhill, reduce speed in accordance with the degree of the incline.
- Do not cross differences in height of more than 5 cm.
- Do not cross steps or curbs without braking first.

The power wheelchair has been approved for driving on inclines and slopes of no more than 17 %. You are not allowed to drive on inclines or slopes steeper than this.

To ensure safe downhill driving, the driving speed must be reduced in accordance with the degree of incline (e. g. set it to speed level 1).

The critical obstacle height of the power wheelchair is 5 cm. Differences in height greater than 5 cm must not be crossed. Do not cross obstacles such as steps or curbs without braking first.

If there are obstacles in the wheelchair's way, drive around them. On uneven ground, driving behaviour of the wheelchair may get out of control. For this reason, the speed of the power wheelchair must always be adapted to the ground conditions.

In the driving menu, the battery status and speed level are indicated on the display.

The joystick (see fig. 13, item 1) is used for driving. The further the joystick is moved away from the mid-position, the faster the power wheelchair will drive in this direction.

The maximum speed with full deflection of the joystick is dependent on the selected speed level.

Letting go of the joystick automatically activates the brake function, which brings the wheelchair to a halt. When standing still, the mechanical brakes are automatically active so the power wheelchair cannot move.

The power wheelchair has 5 speed levels. The button "M" (see fig. 13, item 6) is used to increase the speed level. After speed level 5, the control unit switches back to speed level 1.

Adapting the Driving Characteristics

DANGER!

Risk of accident and injury as a result of incorrect programming!

Programming bears the risk of faults, which can result in uncontrolled driving behaviour. For this reason, test the driving behaviour of the power wheelchair after programming.

A hand programming device, which is to be connected to the control panel, can be used by service staff to adapt the speed, acceleration and deceleration values to the individual wishes of the user.

7.3.4 Battery Indicator

The battery indicator on the LCD display consists of 7 segments and shows the remaining battery capacity.

NOTE!

Immediately after switching the power wheelchair on, the battery indicator shows the battery capacity saved before the wheelchair was switched off the last time. The exact state of the remaining capacity will be displayed after driving for a short time.

A charge of 100% corresponds to 7 full segments on the battery symbol. If one segment disappears, this means that the battery's capacity has decreased by approximately 14%. The distance range on level ground is approx. 35 km, which means that each segment corresponds to a distance range of approximately 5 km if the speed is kept constant.

If the last segment flashes, the battery must be charged immediately.

If all the segments have disappeared and the battery symbol is flashing, the battery is under voltage. Since any further use will lead to battery damage, a warning signal also appears.

If the battery symbol and all the segments are flashing, this means that the battery is over voltage.

Since this will lead to battery damage, a warning symbol also appears.

The charging process is indicated by the battery segments lighting up one after the other. When the battery is charging, the driving function is blocked.

Display symbol(s)	Information
	Driving menu with speed level and battery capacity
	Low battery capacity
	Charging process with drive stop
	Battery is under voltage, with warning symbol
	Battery is over voltage, with warning symbol

Table 3 Battery indicator on the control panel

7.3.5 Electric Seat Functions

NOTE!

Please note that your power wheelchair may always be equipped with a **maximum of two** electric seat functions at a time. Should the electric footrests be controlled individually, they operate as two independent electric seat functions. If the electric footrests, left/right can be controlled jointly, they operate as an electric seat function.

The electric seat functions, such as electric back angle adjustment or seat tilt, are also activated by pressing the button “M” (approx. 2 seconds). Moving the joystick to the right lets you switch between the different functions. The function that is currently selected is shown on the LCD monitor.

The respective function can be extended or retracted by moving the joystick ahead or back. The electric motor moves the respective function as long as the joystick is deflected, and stops at the end positions.

Electric back angle adjustment: When the joystick is deflected ahead or back, the back rest moves ahead or back accordingly.

Electric seat tilt: Deflecting the joystick to the rear will electrically tilt the seat to the rear. Deflecting the joystick to the front tilts the seat ahead until it is level (factory setting; can be adapted by the specialist dealer).

Lift seat (only for B500-S available): Deflecting the joystick ahead or back raises or lowers the electric lift seat. As soon as the lift seat leaves the lowest position, the speed is reduced (creep speed) because the power wheelchair is less tip-resistant. This is indicated by a snail symbol on the control panel LCD monitor. (see Figure 24).

Electric footrests: Deflecting the joystick ahead or back raises or lowers the electric footrests.

Coupled seat adjustment: Deflecting the joystick ahead or back moves the back rest and seat forward and down or back and up simultaneously.

You can switch back to driving mode by pressing the “M” button.

Display symbol(s)	Information
	Electric back angle adjustment
	Lift seat (with B500-S only)
	Electric seat tilt If the lift seat and electric footrests are ordered together, this symbol is used for the lift seat.
	Electric footrests
	Coupled electric back angle adjustment and seat tilt

Table 4 Display of electric seat functions on the control panel

7.3.6 Drive-away Lock

The control unit of the B500 power wheelchair features an electronic drive-away lock for blocking the driving function.

Activating the drive-away lock:

1. If the control unit is turned on, press and hold the On/Off button for at least 2 seconds.
2. A short beep confirms the activation of the drive-away lock (if you switch on the feature ensuring a beep with each push of the button, then it is the second beep that confirms the activation of the drive-away lock).
3. The control unit is turned off automatically. Activation of the drive-away lock is indicated by the key symbol on the display.

Display symbol(s)	Information
	Drive-away lock

Table 5 Display of the drive-away lock on the control panel

Deactivating the drive-away lock

1. Press the On/Off button to turn on the control.

2. Move the joystick into the front end position until you hear a beep.
 3. Move the joystick into the rear end position until you hear a beep.
 4. Release the joystick.
 5. Another beep confirms the successful deactivation of the drive-away lock.
 6. Speed level and battery indicator appear on the LCD display, and the control unit is ready for driving.
- Now the drive-away lock is deactivated and it is possible to drive.

NOTE!

If the joystick is not moved correctly, the drive-away lock remains activated. Turn the controls of the power wheelchair off and on again and make another attempt at deactivating the drive-away lock.

7.4 Releasing and Locking the Brake

If the control unit fails or there is not enough battery power, it is possible to push the power wheelchair by disengaging the brake by means of the mechanical release mechanism. The brake release mechanism is located either to the right or left side of the wheelchair between the frame and seat.

**DANGER!**

Danger to life if the brake fails!
An incorrectly set brake can lead to brake failure and thus result in severe injuries or death. Repairs and adjustments to the brake may only be carried out by authorised service staff.

**DANGER!**

Risk of accident and injury if brake function is unavailable!
When the brake is released (push mode), no brake function is available. When moving on inclined paths, the person pushing the power wheelchair must provide the appropriate brake force.

**DANGER!**

Releasing the brake may result in uncontrolled rolling of the power wheelchair. Therefore make sure that the brake is engaged after parking the power wheelchair.

Releasing the Brake (Fig. 15)

For reasons of safety you must first pull up the brake release bolt.

**NOTE!**

When in this position, the brake must still be engaged! Should the brake already be released and the pushing function be active at this stage, please contact your wheelchair supplier.

Now the brake release lever has been released and can be pressed fully forwards and downwards.

In this position, the control unit will recognise that the brake has been released and automatically deactivate the driving function. As soon as the joystick is deflected, a warning will appear on the LCD display.

**NOTE!**

After turning the brake release lever, all brake systems are deactivated.

Fig. 15 Brake release

- 1 Brake release bolt
- 2 Brake release lever

Locking the Brake

To lock the brake, push the brake release lever upwards until the brake release bolt engages (see fig. 16). To reactivate the driving function, switch off the wheelchair controls and then switch them on again.

Fig. 16 Locking the brake

7.5 Batteries

The standard version of the B500 power wheelchair includes two low-maintenance 12 V wet cell batteries.

NOTE!

The liquid level in the batteries must be checked on a monthly basis and be topped up with distilled water if required.

The batteries are located in a battery case under the seat of the power wheelchair.

To check the batteries, proceed as follows:

**NOTE!**

When checking the power wheelchair's batteries, make sure the power wheelchair is standing on level ground.

**NOTE!**

Before doing any maintenance work on the batteries, please read the battery manufacturer's warnings thoroughly!

1. The seat bottom of the wheelchair must be removed for maintenance work and for mounting/removing the batteries.
For this purpose, undo the hook and loop connection between the backrest and seat bottom and remove the seat bottom.
2. Undo the snap lock on the safety strap for the battery case (see fig. 17) and remove the cover. Once the cover has been removed, you have free access to the batteries (see fig. 18).

Fig. 17 Snap lock on safety strap

3. Unscrew the caps (see arrow in fig. 18) of the individual chambers by using a large screwdriver. You will be able to see the battery's liquid level on the acid level markings.
4. If the liquid is below this marking, top up with distilled water.


Fig. 18 Batteries

5. Close the battery case and reinsert and attach the seat bottom.

Maintenance-free lead gel batteries can be used as an alternative.

See section 7.3.4 for battery capacity display on the control panel.

7.5.1 Charging

DANGER!

Explosion hazard!

Explosive gases can develop while the batteries are charging. The following safety instructions must be followed under all circumstances:

Ensure sufficient ventilation when charging the batteries in closed rooms.

Smoking and fire are not permitted.

Sparks must be avoided.

Do not cover the vent slots in the housing.

ATTENTION!

Battery damage!

Driving over a longer period of time with only the lower segments of the battery indicator lit will discharge the battery completely and damage it. There is a risk that the power wheelchair may stop due to zero battery capacity and bring the user into a dangerous situation.

The remaining battery capacity determines the distance range of the power wheelchair. Many factors have an influence on the battery capacity. Besides the temperature, the age of the batteries and the time and extent of

driving, the way you charge the batteries also has a considerable effect on their capacity and thus on the distance range. For an optimal charging frequency the following applies:

- The batteries can be charged at any time independently of their remaining capacity.
- It takes about 10 hours until a discharged battery (only one, flashing segment) is completely charged. When the charging process is complete, the battery charger can remain connected with no risk of overcharging or damaging the battery. The battery charger features a programmed recharging phase that will maintain the battery capacity at the previously reached level.
- If the power wheelchair is used daily, we recommend connecting the battery charger over night to have the full capacity available every day.
- When the power wheelchair is not used for an extended period of time, the battery will become discharged over time. It is therefore necessary to charge the batteries at least once a week to maintain full battery capacity. We also recommend removing the fuse, when the power wheelchair is not going to be used for an extended period of time.
- Total discharge of the batteries should be avoided.
- When the batteries are charging, the controls on the power wheelchair must be switched off to allow the

charging current to be fed completely into the battery. For battery charging the following must be observed:

- Only the battery charger provided by Otto Bock may be used for charging. Failure to comply will render the warranty null and void.
- The voltage settings on the battery charger must correspond to the voltage used in your country.

7.5.2 Battery Charger

The battery charger is designed for maintenance-free and low-maintenance batteries. The charger can operate along two characteristic curves to charge the respective batteries. The power wheelchair is delivered with the appropriate characteristic curve, set in accordance with the kind of battery the wheelchair is equipped with. If the battery charger is to be used with another power wheelchair or if new batteries are installed, the characteristic curve must be checked.

ATTENTION!

Battery damage!

Setting the characteristic curve incorrectly can result in permanent damage to the battery.

**ATTENTION!**

Unauthorized battery replacement!
The battery may only be replaced by a specialist dealer. The characteristic curve of the battery charger set at the factory corresponds to the battery provided and must not be changed. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

The battery type can be set by means of a switch located at the bottom of the battery charger (covered by a green adhesive label). Once the battery type has been set, the green LED flashes at the beginning of the charging process.

When handling the battery charger, the following safety instructions must be observed:

- Always place the rubber feet of the battery charger on a level surface.
- Protect the battery charger from direct sunlight to prevent the charger from getting any warmer than it already is.
- The place where the battery charger will be used must be dry and well ventilated. Avoid letting dust and dirt enter and affect the function of the battery charger.
- To clean the charger, use a dry piece of cloth.

**DANGER!**

Explosion hazard as a result of sparks!
Always switch off the battery charger and remove the mains plug before disconnecting the battery.

To charge the batteries, proceed as follows:

1. Turn off the controls on the power wheelchair.
2. Connect the plug of the battery charger to the charging receptacle on the control panel of the power wheelchair (see section 7.3.1, fig. 13, Pos.10).
3. Connect the battery charger to a wall socket and turn it on. The batteries will start charging automatically, and the current charge status is indicated on the LEDs of the battery charger (see table 6).
4. When the charging process is complete, turn off the battery charger: Disconnect the battery charger from the mains supply, then from the control panel.
5. Turn on the controls, and the power wheelchair is ready for driving.

The following states are indicated by the battery charger:

Indicator	Function
Yellow LED is lit	Battery is charging
Green LED flashes 1 x	Low-maintenance battery
2 x	Maintenance-free battery
if permanently lit	Battery is charged to 80%
Green LED is lit	Battery is fully charged
Red LED is lit	Defective battery; battery is missing; charging time has been exceeded

Table 6 Battery charging states during the charging process

If none of the LEDs are lit, there is no power supply.

8 Options

ATTENTION!

Only original options provided by the manufacturer may be used. The optional components may be mounted only as described here. Failure to comply will render the warranty null and void.

NOTE!

All available optional components are contained on the order form and in the wheelchair accessories catalogue.

The B500 power wheelchair has been designed as a modular system. Certain component groups can be exchanged and other accessory components can be added to the power wheelchair. The complete range of accessory components is listed on the order form and in the wheelchair accessories catalogue. Some important options are described in the following.

8.1 Seat Options

DANGER!

Risk of injury due to tipping during operation. Before driving up slopes or over obstacles, a tilted seat must be lowered to its basic position, the lift seat must be in its lowest position, and the back must be in the upright position. Slightly tilt the seat to the rear when driving down inclines or over obstacles and reduce the speed (max. 3 km/h).

DANGER!

Risk of injury if the power wheelchair tips over! Before driving up slopes or over obstacles, a tilted seat must be lowered to its basic position and the back be brought into the upright position.

DANGER!

Risk of pinching or crushing the fingers! Do not reach into the danger zone between the frame and the seat when tilting the seat. There must not be any interfering objects or obstacles in the tilting area.

DANGER!

Overloading the actuators can result in spindle nut breakage and thus cause the seat to drop or the backrest to fold back.

NOTE!

It is easier to tilt the seat, if the backrest is inclined to the rear.

NOTE!

The actuators of the seat functions are not designed for continuous operation but rather for an operational load that is limited to a short time only (10% operation - 90% stop). General recommendation: With maximum load, the actuators may be activated for 10 seconds; then they should not be used for approximately 90 seconds. The electric seat functions are not considered part of the driving function in this case.

NOTE!

When getting into or out of the power wheelchair, bring the seat to a level position.

8.1.1 Mechanical Seat Tilt Adjustment

The mechanical seat tilt adjustment lets you bring the seat into a tilted position for pressure relief, for example. The standard or contoured seat can be continuously tilted backwards up to 20°.

Fig. 19 Mechanical seat tilt adjustment

DANGER!

Risk of injury as a result of abrupt lowering!
Using the release lever when the seat is tilted can lower the seat abruptly.

To tilt the seat backwards (see fig. 19):

- Use the release lever on the armrest (to unblock the gas compression spring).

- Tilt the seat to the desired inclination.
- Let go of the release lever (to again block the gas compression spring).

8.1.2 Electric Seat Tilt Adjustment

The electric seat tilt adjustment lets you bring the seat into a tilted position of up to 22° for pressure relief, for example. To control the tilting, the joystick is used in the “additional function” mode. The seat can be continuously tilted backwards (see fig. 20).

NOTE!

It is helpful to slightly tilt the seat to the rear when driving down inclines or curbs. Drive down inclines and across obstacles at reduced speed only (max. 3 km/h).

Fig. 20 Electric seat tilt adjustment

8.1.3 Mechanical Back Angle Adjustment

The backrest of the B500 power wheelchair with a standard or contoured seat can be equipped with a mechanical back angle adjustment.

To adjust the backrest (see fig. 21):

- Use the release lever on the armrest (to unblock the gas compression spring).
- Move the backrest to the desired position.
- Let go of the release lever (to again block the gas compression spring).

Fig. 21 Mechanical back angle adjustment

8.1.4 Electric Back Angle Adjustment

The backrest of the standard seat can be equipped with an electric back angle adjustment (see fig. 22). To control this feature, the joystick is used in the “additional function” mode. The backrest can be continuously inclined backwards up to 30°.

NOTE!

It is helpful to slightly incline the back to the rear when driving down inclines or curbs. Drive up and down inclines, slopes, and obstacles at reduced speed only (max.3 km/h).

Fig. 22 Electric back angle adjustment

8.1.5 Lift Seat Function (B500-S only)

DANGER!

Risk of injury if the power wheelchair tips over. The power wheelchair can tip over when the lift seat is raised, causing the user to fall out of the wheelchair.

The following safety instructions must be followed under all circumstances:

- Only travel short distances when the lift seat is raised, e. g. when manoeuvring.

- Use the lift seat function only with the back rest set to a vertical position.
- Apply the lap belt if possible.
- Only use the lift seat function on level surfaces!
- Do not change the mounting position of the lift seat mechanism or the lift seat.
- For lift seat retrofit or conversion, observe the instructions provided in the service manual.

DANGER!

Risk of injury if the power wheelchair tips over. Creep speed is automatically activated when the lift seat is operated. If this is not case, it represents a malfunction and the lift seat function must not be used. Contact a specialist dealer as soon as possible to eliminate the problem.

DANGER!

Risk of pinching in the lifting zone. When using the lift seat function, no unauthorised persons may be present in the danger zone. There must not be any interfering objects or obstacles in

the lifting zone.

All attendants must be informed that there are pinch points caused by the construction in the area between the seat frame and wheelchair frame. Users and attendants must not reach into the danger zone if the electric seat tilt and / or lift seat functions are installed.

**DANGER!**

Risk of falling. Do not lean the upper body beyond the edge of the seat when the lift seat is raised.

**DANGER!**

Risk of accidents due to a restricted field of vision. Since the user's field of vision is restricted when the seat is raised, extra caution is required when operating the wheelchair with the lift seat in the raised position and when activating the lift seat function.

**DANGER!**

Risk of tipping in case of additional use of the seat tilt function. The seat tilt function must be set to the horizontal position when driving

up inclines or over curbs. Since tipping resistance is reduced in such situations, additional caution is required!

**ATTENTION!**

Risk of damage due to improper maintenance. The lift seat must be checked for visible damage and the tightness of all screw connections must be verified at least once a month. Also maintain correct air pressure.

**ATTENTION!**

Risk of damage during transportation. Always lower the lift seat fully for loading or transportation.

**ATTENTION!**

Risk of damage due to overloading. The maximum load capacity of the lift seat is 130 kg.

**NOTE!**

A sufficient supply of power must be maintained in order to guarantee proper operation of the lift seat.

**NOTE!**

The lift seat is not intended for continuous use, but only for short-term limited operation (10 % use, 90 % idle time). For every 1 hour of operation, the lift seat function can be used for a maximum of 6 minutes.

The B500-S can be equipped with an optional lift seat. With the lift seat function, the seat can be raised by up to 40 cm using a motor (Figure 23). It is controlled by the joystick in “additional function” mode.

The drive function remains available even when the lift seat is raised. As soon as the lift seat leaves the lowest position, the speed is reduced (creep speed) because the power wheelchair is less tip-resistant. This is indicated by a snail symbol on the control panel or the optional, separate LCD monitor (Figure 24).

Fig. 23 Lift seat in raised position (B500-S)

Fig. 24 Snail symbol indicating creep speed

8.1.6 Preparing Wheelchairs with Back Angle Adjustment for Transport

To obtain the folding size, the following must be done with the back angle adjustment:

- Release the cross bolt at the lower end of the gas compression spring (by pushing up the lever, see arrow in fig. 25).
- Remove the bolt from the holding device.
- Remove the side panels and fold down the back part onto the seat bottom.


Fig. 25 Cross bolt locked in the holding device

- 1 Cross bolt
- 2 Lock with lever

To put a wheelchair with back angle adjustment into operation:

- Unfold the backrest.
- Insert the cross bolt at the end of the gas compression spring into the holding device.
- Push down the lever to lock the bolt (see fig. 25).

ATTENTION!

Verify that the bolt and lock of the back angle adjustment are secure!

8.2 Lighting

An option for the power wheelchair is to have it equipped with a set of lights.

The set of lights consists of:

- Two rear lamps with integrated direction indicators
- Two front lights consisting of halogen lamp and direction indicator (fig. 26).

Fig. 26 Front lights on the side panel

The front lights are mounted on the side panel. The rear lights are attached to the wheelchair frame.

Button 9 on the control panel (see fig. 13) is used to switch the lights on.

Symbol 16 on the LCD display (see fig. 14) indicates that the lights are turned on.

Changing Defective Lights

NOTE!

Lamp housings or lamps can be ordered from a specialist dealer.

To change the halogen lamp in the front light, proceed as follows:

- Lightly push the black engaging lever at the lower end of the lamp and fold the pane forward (see fig. 27, item 1).


Fig. 27 Changing the front lights

- 1 Engaging lever of front light
- 2 Recess (use the screwdriver at this point)

- Remove the pane with mirror towards the front.
- Pull the bulb out of the rear of the mirror.
- Replace and reinsert the bulb.
- Insert the pane with mirror and bulb into the lamp housing.
- Make sure that the pane securely snaps in place.

To change the bulb of the front direction indicator proceed as follows:

- Open the front direction indicator: Insert a small screwdriver into the recess of the direction indicator housing and flip down the pane of the direction indicator (see fig. 27, item 2).
- Remove the pane.
- Slightly rotate the bulb to release it from the lock on the lamp socket and pull it out (see fig. 28).

Fig. 28 Replacing the direction indicator lamp

To change the rear direction indicators and the rear light, proceed as follows (see fig. 29):

- Unscrew the pane.
- The rear lamps are attached on the right and left side on sprung-mountings. Exert slight pressure to the right and pull the lamp forward on the left side.

Fig. 29 Rear light without lamp cover

To protect the front lamps against damage, they are provided with rubber bearings.

ATTENTION!

Do not expose the lamps to moisture!
When reinserting the panes, make sure they sit correctly on the housing and that the screws are firmly tightened.

8.3 Suspension of Casters and Drive Wheels

The power wheelchair can be equipped with a suspension system. One spring/damping element on each of the two casters increases driving comfort and provides better adherence to the ground, especially when driving on uneven terrain. This can be improved even more by also using spring-mounted drive wheels.

To order a suspension system, a specialist dealer needs to know the user's weight in order to determine the degree of spring damping.

When equipped with drive wheel suspension, the B 500 has a battery drawer.

To open the battery drawer:

- Pull up the release bolt (see fig. 30).
- Pull the drawer out backwards.
- Undo the snap lock on the safety strap.
- Remove the cover of the battery case.

Fig. 30 Release bolt on the battery drawer

8.4 Front Wheel Swivel Lock

The wheelchair must be driven in totally straight line onto ramps or into lifts. The B500 power wheelchair also offers the option of having the front wheels locked in the forward position, thereby preventing the wheelchair from turning.

Engaging the front wheel swivel lock:

- Slightly rotate the folding lever on the front frame to the side thus allowing it to jump out of its normal mid-position.

- A bolt is released from the frame and engages in the front caster fork as soon as the straight forward position of the wheelchair is reached (see fig. 31, left photo).

Now the power wheelchair can drive straight forward or backward.

Fig. 31 Caster swivel lock engaged and released

Releasing the front wheel swivel lock:

- Turn the folding lever back to its mid-position where it will disengage the bolt from the front caster fork.
- The lever engages into the centre of the front frame in an unlocked position (see fig. 31, right photo).

8.5 Elevating Footrests

8.5.1 Mechanically Elevating Footrest

The B500 power wheelchair can be equipped with mechanically elevating footrests (see fig. 32).

To elevate the footrest:

- Use the release lever on the footrest (to unblock the gas compression spring).
- Move the footrest to the desired position.
- Let go of the release lever (to again block the gas compression spring).

Fig. 32 Mechanically elevating footrest

8.5.2 Electrically Adjustable Footrests

NOTE!

The seat function actuators are not intended for continuous use, but only for short-term limited operation. General recommendation: At maximum load capacity, 10 seconds of operation should be followed by 90 seconds of idle time. The electric seat functions are considered independently of the driving function for this purpose.

To avoid permanent pressure loads or to ensure anti-shock positioning, the power wheelchair offers the possibility of mounting electrically adjustable footrests (Figure 33). In the control unit's „additional function“ mode, both footrests can be controlled separately or simultaneously depending on the ordered version. To adjust them, use the joystick.

The amount of space available for getting into or out of the wheelchair can be increased by flipping the footplates up. (Figure 11).

Fig. 33 Elevated footrest

8.6 Control Panel Holders

8.6.1 Height Adjustable Control Panel Holder

Another option is a special control panel holder that can be mounted on the B500 power wheelchair, allowing the control panel to be moved downwards. To make use of this option, loosen the attachment screw (see fig. 34) and move the control panel to the desired height. Firmly retighten the attachment screw.

Fig. 34 Height adjustment of control panel holder

8.6.2 Swing-away Control Panel Holder

To allow the user to drive the power wheelchair closer to an object or under the edge of a table, the control panel can be swung to the side with a special control panel holder (see fig. 35).

The control panel holder can be turned up to the armrest:

- Apply a little bit of pressure to press the control panel holder to the side and release the pivot element.
- Swing the control panel holder to the side.
- When bringing the control panel holder back to its original position, the pivot element automatically locks into place.

Fig. 35 Swing-away control panel holder

8.7 External Mileage Counter

The external mileage counter is attached to the protective bow for the control panel.

The mileage counter shows:

- Speed
- Trip mileage
- Total mileage
- Digital clock

Fig. 36 External mileage counter

Whenever the display shows km/h or m/h in the upper left corner, the speedometer is activated.

Pressing the yellow button lets you see the trip mileage, total mileage, digital clock and speed one after the other on the display.

Setting the Trip Mileage Back to Zero

Three flashing arrows on the left side of the display indicate the trip mileage counter. To set the trip mileage back to zero, press the yellow key for at least two seconds.

Setting the Time

If the digital clock is activated on the display, you may set the time. Press the yellow key for at least two seconds. The display will show 12: and 24: in turn. Pressing the key again when the display shows 24: will program the clock to show a 24-hour display. Pressing the key when the display shows 12: will result in a 12-hour display.

Next, the hours and minutes can be set. The time digits automatically advance one after the other. Press the yellow key when the desired digits are showed on the display.

Setting the Total Mileage

The mode for the total mileage counter lets you switch from a mileage to kilometre display and enter the wheel circumference. To do so, press the yellow key for at least 2 seconds and then again press the key to select between the alternating displays (km/h and m/h). The wheel circumference can be calculated on the basis of the wheel size and can be set in the same way the clock is.

For a wheel size of 14 inches:

Conversion from inches to cm	14" x 25.4 = 355.6 mm
$C = d \times \pi$	355.6 mm x 3.1416 = 1 117 mm (wheel circumference to be set)

8.8 Recaro® Seat

The B500 power wheelchair can be equipped with different Recaro® seat models.

To adjust the back angle of a Recaro® seat model, use the turning handle on the right or left side of the backrest.

DANGER!

Risk of injury if the power wheelchair tips over! Before driving up slopes or over obstacles, a tilted seat must be lowered to its basic position and the Recaro backrest be brought into the upright position.

DANGER!

When standing still, the backrest may be inclined to the rear up to 30° only. During driving, the backrest may be inclined to the rear up to 20° only.

To fold down the backrest, you can use a release handle that is located laterally at the backrest. By pulling the release handle upward, the backrest can be folded down to the front.

If the backrest is brought into its initial position, it automatically locks in place. Verify that the backrest is securely locked in place every time it has been unlocked!

With the Recaro LT backrest, the lateral guide in the lumbar section can be individually adjusted by turning the hand wheel. Turning the wheel to the front brings the two lateral guides closer to each other evenly. Turning the wheel to the rear adjusts the two lateral guides away from each other.

Seat Bottom

The seat bottoms X and W have a telescoping padding section at the front. The lever used to adjust the padding section is located below the seat bottom. Pull the lever up to shift the padding section into the desired position. The padding section locks in place as soon as the adjustment lever is released. The adjustment lever must audibly lock in place and return to its initial position.

Headrest

A fixed height of the headrest can be set; the inclination is adjustable.

A second person is needed to remove the headrest. The pressure points must be detected under the cover by feel and be pressed simultaneously. Now the second person can pull the headrest out upwards.

Removing the Recaro® Seat from the Chassis

- Pull the release strap beneath the front part of the seat.
- Loosen the right and left bolt that serve to lock the seat to the frame.
- To release the posterior seat attachment from the retainer bushings, slightly tilt the seat to the rear and slide back a little.
- Remove the Recaro® seat.

Attaching the Recaro® Seat to the Chassis

ATTENTION!

Risk of injury if lock does not engage!
When mounting the Recaro® seat, the seat adapter must be securely attached to the retainer bushings. The anterior locking bolts must have snapped into place up to the key ring.

- Place the posterior edge of the Recaro® seat on the end of the seat frame.
- Slide the seat forward until the posterior seat attachment hangs in the retainer bushings.
- Tilt the seat to the front until the locking bolts engage with the seat frame.

8.9 Other Options

8.9.1 Lap Belt

A lap belt can be added to the power wheelchair. The lap belt helps to additionally stabilise the person sitting in the wheelchair and prevents him or her from slipping out of the seat.

DANGER!

The lap belt helps to additionally stabilise the person sitting in the power wheelchair. The lap belt must never be used as part of a restraint system for transportation of the person in motor vehicles for the disabled.

Fig. 37 Applying the lap belt

To apply the lap belt, insert the two buckle halves into each other until they lock in place. The belt buckle must snap in place audibly. Then, verify that the belt has locked by trying to pulling it apart. The lap belt should not be too tight on the body. Objects caught under the belt can cause painful pressure sores.

To open the lap belt, press the red release buttons.

The belt length can be adjusted on both sides. To set the length, position both buckle halves in the middle of the body. After positioning the buckle halves at a right angle in relation to the belt, their position can be varied. Any excess belt length will be held in place by the plastic slides.

8.9.2 Four-way Chest Strap Kit

The four-way chest strap kit makes it possible to restrain the patient in the power wheelchair. The lower strap should be placed between the lap belt and thigh, and should not be too tight. The belt buckle is positioned at the front of the body, in the centre. The shoulder straps should be positioned evenly on both sides on the thorax, and should not be too tight.

**DANGER!**

The four-way chest strap kit helps to additionally stabilise the person sitting in the power wheelchair. The four-way chest strap kit must never be used as part of a restraint system for transportation of the person in motor vehicles for the disabled.

To close the four-way chest strap kit, connect all parts of the belt buckle. The belt buckle must snap in place audibly. Then, verify that the belt has locked by trying to pulling it apart. The four-way chest strap kit should not be too tight on the body. Objects caught under the belt can cause painful pressure sores.

To increase the length of the belt, turn the buckle by 90° and pull. To reduce the length of the belt, pull the overhanging belt end.

8.9.3 Headrest Mounting Kit (Standard Seat)

Another option is to have the power wheelchair equipped with a headrest. For this purpose, the mounting kit showed in figure 38 is required.

Fig. 38 Headrest mounting kit

8.9.4 Rear Bumper

Equipping the B500 power wheelchair with a bumper is also an option (fig. 39).

Fig. 39 Rear bumper

8.9.5 Other Options

- Curb climbing assist
The curb climbing assist is used to clear curbs and steps that have a maximum height of 10 cm.
- Puncture-proof tyres
Solid rubber tyres
- Contoured seats are available in different sizes for better sitting comfort and more adaptation possibilities.
- Seating shell interface
For the adaptation of special seating shells. Control panel holders for seating shells are also available.

- Armrest accessories
Special adapter for the armrests in our wheelchair accessories catalogue
 - Control for attendant
External joystick on the backrest
 - Joystick accessories
Fork for tetraplegics, STICK S80, soft ball, golf ball, flexible control stick shaft
 - Protective bow for control panel
Metal bow for protecting the control panel against shocks and damage
 - Crutch holder
 - Folding rear-view mirror
 - Tray
Attachable tray top
 - Accessories case, backpack, pocket for mobile phone
- These and other optional add-on components are included in the order form and in the wheelchair accessories catalogue.

9 Failures/Troubleshooting

NOTE!

Should you encounter problems while troubleshooting or if you do not manage to completely eliminate a problem by following the measures described here, please contact your specialist dealer.

Problems are indicated on the display of the control panel. Table 7 describes the individual indications with the corresponding sources of the problem as well as the possible causes and measures to be taken.

If you do not manage to completely eliminate the problems with the described measures, your specialist dealer may be able to read the exact error code with the hand programming device and undertake a more systematic analysis of the system.

All problems that have ever occurred are saved in a list and can be retrieved, for example, in case of a general overhaul of the power wheelchair. The saved data can be used for determining future service and maintenance intervals, for example.

9.1 Warning

A warning indicates a status or malfunction of one or several individual components of the power wheelchair. The function of components without errors is not restricted. If, for example, the connection between the controller and seat motor malfunctions, this error will only be indicated if the motor is used. The driving function is still available, however.

9.2 Error

ATTENTION!

Risk of injury if the power wheelchair stops abruptly!

In the event of communication problems in the bus system of the controls, the system triggers an emergency stop to prevent any uncontrolled functions. After switching the wheelchair controls on again, depending on the problem, the wheelchair can be driven out of a danger zone, e.g. a crosswalk. If the driving function is still not available after switching on the controls again, the brakes must be released (see section 7.4) to switch over to the pushing mode. After the power wheelchair has triggered an emergency stop, contact a specialist dealer as soon as possible!

ATTENTION!

Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunction. In this case, please contact your authorised dealer immediately.

An error affects one or several functions of the power wheelchair. The system is not fully operational until the error has been eliminated.

9.3 Defect/Failure

A defect represents a serious failure of a system component.

After elimination of the fault, the system will be re-activated by re-starting it.

A failure represents the most serious kind of fault possible and results in an immediate emergency stop of the system.

Display symbol(s)	Error/warning	Cause
	Communication error (flash alternatingly)	Faulty connection between hand control device and controller/ defective cabling, software or hardware
	Battery: under voltage	Total discharge of battery/charge as soon as possible
	Battery: over voltage	Voltage too high (after completed charging process and downhill drive)
	Fault with seat tilt motor	Faulty cabling or plug contact, defective actuator

Display symbol(s)	Error/warning	Cause
	Controller temperature warning	Overheating due to excessive load >> cooling down phase
	Motor temperature warning	Overheating due to excessive load >> cooling down phase
	Joystick warning	Joystick not in zero position when switching on
	Fault with hand control device	Defective joystick
	Controller fault	Defective controller
	Fault with back adjustment motor	Faulty cabling or plug contact, defective actuator

Display symbol(s)	Error/warning	Cause
	Drive motor fault	Faulty cabling or plug contact, defective actuator
	Brake fault	Open brake release/defective brake
	Emergency stop	Severe fault caused by malfunction in the controller, hand control device and/or drive motor

Table 7 Status and error messages

10 Maintenance, Cleaning and Care

**NOTE!**

A spare parts catalogue from Otto Bock is available upon request for ordering spare parts. Only spare parts and replacement parts from Otto Bock may be used. Failure to comply will render the warranty null and void.

**NOTE!**

Should you encounter problems during maintenance, contact your authorised dealer. The power wheelchair is to be checked and serviced once a year by an authorised dealer. Tyre changes and repairs are also carried out by the specialist dealer.

10.1 Maintenance Intervals

The correct function of the power wheelchair should be checked every time before using it. The items listed in table 8 must be checked by the user at the indicated intervals.

Component	Function or part	Daily	Weekly	Monthly
Drive wheels	Wheels must rotate freely and without axial run out			X
	Central nut on the driving shaft tightened			X
	Check if wheel mounts are seated securely			X
	Directional stability of the entire wheelchair		Prior to every ride	
Steering wheels and casters	Fork seated in the receiver without play			X
	Wheels must rotate freely and without axial run out			X
	Fastening nut tightened			X
Seat attachment	Check if mounting screws are securely in place			X
	Check seat lock of Recaro® seat			X
Footrest	Check ratchet mechanism for functionality and secure seating			X
	Check footrest for damage			X
Mechanically elevating footrest	Check ratchet mechanism for functionality and secure seating			X
	Check footplates for damage			X
	Visual check for scratches on the piston rod and oil leakage			X
Padding and belts	Proper condition of padding			X
	No wear on the seat belts			X
	Check belt buckle for functionality		X	
Tyres	Air pressure (printed on the sidewall of the tyre)			X

Component	Function or part	Daily	Weekly	Monthly
Tyres	Sufficient tread depth, at least 1 mm Check for damage			X
Batteries	Check liquid / acid level (not with gel batteries)			X
Lighting	Perform visual check for damages Test to see that it functions	X	X	
Electronics	Control without error message Battery charger does not show any error messages on the LEDs Check plug connections		Prior to every ride X	 X
Brake	Use the brake lever with the controls turned on Brake function active with engaged brake	X		 X
Lift seat	Visual inspection of all moving components, especially cabling - check for damage Check whether screw connections are tight			X X
Caster stabilisation	Clean and grease the entire component group			X
Side panel and armrest	Fastening screw tightened Armrest and control panel secured Check armrest for damage	X	X	X
Gas compression spring or actuator	Visual examination for scratches on the piston rod and oil leakage			X

Table 8 Maintenance measures and intervals

NOTE!

If the brake has been released and the joystick activated, the controls will emit an error signal on the control panel. If this is not the case, something has malfunctioned and the problem must be eliminated immediately by a specialist dealer.

10.2 Changing the Fuse

The 80 A fuse is located in a fuse holder at the rear end of the battery case (see fig. 40).

Open the cap of the safety holder and exchange the fuse. Make sure that the fuse is pressed into the centre of the spring contacts and that it is not at an angle. Close the cap until you feel it snap back into place.


Fig. 40 Fuse holder

- 1 Cap, open
- 2 Fuse, inserted
- 3 Fuse holder

10.3 Changing Tyres

Proceed as follows:

- ❑ To dismount a drive wheel, loosen the four Allen head screws in the middle of the wheel (fig. 41) and remove the wheel.

Fig. 41 Disassembling drive wheel

- ❑ If the inner tube of a drive wheel needs to be replaced, loosen all Allen head screws on the inner side of the rim and pull apart the two-piece rim.

The defective inner tube is now freely accessible and can be replaced.

- ❑ To dismount the front wheel, loosen the screw of the axle using a 6 mm Allen wrench (see fig. 42, item 1) and pull out the axle.
- ❑ Loosen the five Allen head screws (see fig. 42, item 2) and separate the two-piece rim.

The defective inner tube is now freely accessible and can be replaced.


Fig. 42 Disassembling front wheel

- 1 Axle screw
- 2 Allen head screw

10.4 Cleaning and Care

ATTENTION!

Risk of damage to the electronics!

Risk of malfunctions!

When cleaning the power wheelchair, the electronic components, the motor and the batteries must not come into contact with water.

ATTENTION!

Risk of damage to components of the power wheelchair!

Do not use any aggressive cleansing agents, solvents, or hard brushes for cleaning the power wheelchair. Never use a water jet or high-pressure cleaning apparatus for cleaning the power wheelchair.

The power wheelchair must be cleaned regularly, depending on the amount of use and the degree of soiling.

Use a damp piece of cloth and a mild cleansing solution to clean the control panel, battery charger, armrest and coverings.

Use a dry brush to clean the seat and back upholstery as well as the seat cushion.

Use a damp plastic brush to clean the wheels and frame.

NOTE!

Prior to disinfection, clean the seat and back upholstery as well as the seat cushion, the control panel and the armrest.

Use a standard disinfection agent like Sagrotan for disinfection.

11 Disposal

**NOTE!**

If a wheelchair is to be disposed of, all components and materials of the power wheelchair must be recycled or disposed of properly.

**NOTE!**

There is a risk of polluting the environment with battery acid!
The batteries of the power wheelchair contain a toxic acid. They must not be disposed of like regular domestic waste and the battery acid must not get into the sewage water system or earth. You must observe the instructions from the batteries' manufacturer, which are printed on the batteries.

If the power wheelchair is no longer in use, it must be disposed of properly in accordance with national regulations.

Please return defective batteries to your dealer when buying new ones.

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Otto Bock HealthCare Canada Ltd.
CA-Canada, Burlington,
Ontario L7L 5N5
Tel. +1 289 288-4848
Fax +1 289 288-4837
info@ottobock.ca

Otto Bock HealthCare Andina Ltda.
Bogotá / Colombia
Tel. +57 1 8619988
Fax +57 1 8619977
ottobock@telesat.com.co

Otto Bock de Mexico S.A. de C.V.
MEX-Guadalajara, Jal. 44210
Tel. +52 33 38246787
Fax +52 33 38531935
info@ottobock.com.mx

Otto Bock HealthCare LP
U.S.A.-Minneapolis,
Minnesota 55447
Phone +1 800 328 4058
Fax +1 800 962 2549
usa.customerservice@ottobock.com

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Otto Bock Mobility Solutions GmbH has been certified by the German Society for the Certification of Quality Assurance Systems (DQS) in accordance with DIN EN ISO 13 485 : 2007, reg. no. 435514 MP27 (management system)

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Otto Bock Korea HealthCare Inc.
Otto Bock South East Asia Co. Ltd.
Otto Bock HealthCare GmbH
Otto Bock Manufacturing Königsee GmbH